KWANTALOG

2017/2018

electronic

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pH - mV - Conductivity - TDS -
Salinity - Temperature

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Application profile

- The ideal handheld instrument to measure pH.
- Quality monitoring in laboratories, productions or stocks
 - pH-measurement in food (meat, sausage and cheese).
 - pH-measurement in soil, aquariums, breweries, laboratories, schools and universities.

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Product features

- Large LCD display
- Buffers used for calibration indication (L= Low, M= Medium, H= High)
- Measurement stability indicator (Smiley-Symbol ☺)
- Automatic temperature compensation
- 1...2 calibration points (PH 1 / COND1) -
1-, 2- or 3-calibration points (PH 5 / PH 5 Food / PC 5)
- Extended conductivity scale: 0...2000µS; 0...20mS
- Conductivity, Salinity and TDS function (COND 1)
- 4x1,5V AAA batteries for long working time, compatible with rechargeable batteries

Detail

Großes Beleuchtetes Farb LCD-Display
Large backlight color LCD display

austauschbarer Se
Replaceable sensc

während der Messung / During measure

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Areas of application

- Agriculture
- Environmental monitoring
- Water and waste water treatment
- Hydroponics Aquaculture and aquariums
- Pools and spas
- Food and beverage manufacturing
- Electroplating
- Printing Education

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5040-0242

Product features

- pH scale: -2...16 pH, Resolution 0.01 pH
- mV function for pH electrode quality test or connection to ORP sensor for Redox potential measurement

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5040-0243

Product features

- with replaceable spear pH electrode. Double open pore polymer electrode maintenance free incorporating temperature sensors for fast temperature reading
- pH scale: -2...16 pH, Resolution 0.01 pH
- mV function for pH electrode quality test

PC 5 pH/mV/Cond/TDS/Salzgehalt/°C / pH/mV/Cond/TDS/Salinity/°C**Produktmerkmale**

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5040-0245

Product features

- pH scale: -2...16 pH, Resolution 0.01 pH
- mV function for pH electrode quality test or connection to ORP sensor for Redox potential measurement
- COND scale: 0...200µS - 0...2000µS - 0...20mS
- Reference temperature 20/25 °C user selectable

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5040-0261

NEU / NEW

DHS Digital High Sensor
Intelligentes Sensor Management
Intelligent Sensors Management

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Product features

- Displays pH / mV / COND and temperature simultaneously
- Changeable resolution 0,1 / 0,01 pH
- Automatic calibration function pH up to 3 calibration points
- Manual calibration pH up to 2 points (user defined buffer)
- Automatic calibration function COND up to 4 calibration points
- Manual calibration COND 1 point (user defined buffer)
- TDS measurement with adjustable factor
- Large lighted display with calibration indicator
- Data logger (automatic or manual) 500 data sets with date & time
- GLP
- CAL DUE: set of calibration frequency for high quality of measurement
- Data Logger, manual and automatic with 500 memory values with date and time
- USB port for data transfer and power by PC or power-supply
- Software PC-Link to download data
- Incl. calibration solution pH 4, pH 7 and storage solution
- IP 57 waterproof
- Dimensions instrument 86 x 196 x 33 mm (L x W x H)
- Weight: approx 300 g
- DHS

Areas of application

- Ideal for the measurement of fluids
- For quality control in laboratories and production process

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Product features

- Displays pH/mV and temperature simultaneously
- Changeable resolution 0,1 / 0,01 pH
- Automatic calibration function up to 3 calibration points
- Large display with calibration indicator
- Incl. calibration solution pH 4, pH 7 and storage solution
- IP 57 waterproof
- DHS

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pH-instrument for measuring of pH, mV and °C incl. special insertion electrode for the measurement of meat, ham, cheese and other foodstuff.

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Areas of application

- In dairy, butchery, brewery and gardening firm

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pH-instrument with data logger for measuring of pH, mV and °C incl. pH electrode 201T with internal temperature probe

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Product features

- Displays pH/mV and temperature simultaneously
- Changeable resolution 0,1 / 0,01 pH
- Automatic calibration function up to 3 calibration points
- Manual calibration up to 2 points (user defined buffer)
- Large lighted display with calibration indicator
- Data logger (automatic or manual) 500 data sets with date & time
- GLP
- CAL DUE: set of calibration frequency for high quality of measurement
- USB port for data transfer and power by PC or power-supply
- Software PC-Link to download data
- Incl. calibration solution pH 4, pH 7 and storage solution
- IP 57 waterproof
- DHS

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Product features

- Displays pH / mV / COND and temperature simultaneously
- Changeable resolution 0,1 / 0,01 pH
- Automatic calibration function pH up to 3 calibration points
- Manual calibration pH up to 2 points (user defined buffer)
- Automatic calibration function COND up to 4 calibration points
- Manual calibration COND 1 point (user defined buffer)
- Large display with calibration indicator
- Data logger (automatic or manual) 500 data sets with date & time
- Incl. calibration solution pH 4, pH 7 and storage solution
- IP 57 waterproof
- DHS

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Product features

- Displays pH / mV / COND and temperature simultaneously
- Changeable resolution 0,1 / 0,01 pH
- Automatic calibration function pH up to 3 calibration points
- Manual calibration pH up to 2 points (user defined buffer)
- Automatic calibration function COND up to 4 calibration points
- Manual calibration COND 1 point (user defined buffer)
- TDS measurement with adjustable factor
- Large lighted display with calibration indicator
- Data logger (automatic or manual) 500 data sets with date & time
- GLP
- CAL DUE: set of calibration frequency for high quality of measurement
- Data Logger, manual and automatic with 500 memory values with date and time
- USB port for data transfer and power by PC or power-supply
- Software PC-Link to download data
- Incl. calibration solution pH 4, pH 7 and storage solution
- IP 57 waterproof
- DHS

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Analysis instrument with logger function and GLP function for measuring of pH, mV/ORP-COND-TDS and temperature incl. Cond cell

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pH
REDOX
CD
TDS
SALT
WIDERSTAND /
RESISTIVITY
-20°C...+120°C
IP 57
BACKLIGHT

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Product features

- Display pH/mV and Temperature, COND/TDS/Salinity/Resistivity and Temperature simultaneously
- Selectable Resolution 0,1 / 0,01 / 0,001 pH
- Auto range for COND
- Calibration for up to 5 pH points with automatic recognition of buffer solution (USA / NIST / DIN)
- manual calibration for up to 5 pH points defined by user
- Calibration for up to 4 COND points with automatic recognition of buffer solution
- Manual calibration for 1 COND point defined by user
- TDS measuring with selectable factor
- GLP and IQ-OQ manual available
- Large color display with indication of Calibration Buffers and measurement's stability
- Memory with re-call of last calibration date and time and information of electrode's life
- CAL DUE: set of calibration frequency for high quality of measurements
- Non-volatile memory holds up to 1000 data points with date and time in automatic (Data logger) or manual and send to PC or printer
- 2 x USB ports to transfer data at PC or Printer and for external keyboard
- Software PC-Link to download data
- Stirrer plate with direct control from PC 80 with independent keys integrated on case or removable with 50cm long cable
- Incl. Buffer solutions and accessories
- IP 54 rated housing handles splashes
- DHS

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Areas of application

- The perfect instrument for Analysis pH/mV, COND, TDS, Salinity, Resistivity and other parameters
- Quality control in laboratories Measurements in water, waste water, HVAC, breweries, dairies and much more products and applications.

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pH-instrument for measuring of pH, mV and °C incl. pH electrode 201T with internal temperature probe and Electrode holder

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Product features

- Display pH/mV and temperature simultaneously
- Resolution 0,01 pH
- Automatic calibration function up to 3 points
- Large display with calibration indication
- Memory function for 25 data points
- IP 54 splash water proof
- Incl. electrode holder
- Incl. pH electrode 201T
- Incl. calibration pH 4, pH 7 and storage solutions
- DHS

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pH-instrument for measuring of pH, mV and °C incl. 1 meter cable for pH electrode S7/BNC, NT55 temperature probe and Electrode holder

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pH bench instrument for measuring pH/ mV and temperature inclusive Electrode holder and accessories

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Product features

- Display pH/mV and Temperature simultaneously
- Selectable Resolution 0,1 / 0,01 pH
- Calibration for up to 3 pH points with automatic recognition of buffer solution (USA / NIST)
- manual calibration for up to 2 pH points defined by user
- Large 2-color display with indication of Calibration Buffers and measurement's stability
- Incl. Buffer solutions and accessories
- IP 54 rated housing handles splashes

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pH bench instrument for measuring pH/ mV and temperature inclusive magnetic stirrer and accessories

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Product features

- Display pH/mV and Temperature, COND/ TDS/Salinity/Resistivity and Temperature simultaneously
- Selectable Resolution 0,1 / 0,01 / 0,001 pH
- Calibration for up to 5 pH points with auto. recognition of buffer solution (USA / NIST / DIN)
- manual calibration for up to 5 pH points defined by user
- GLP and IQ-OQ manual available
- Large color display with indication of Calibration Buffers and measurement's stability
- Memory with re-call of last calibration date and time and information of electrode's life
- CAL DUE: set of calibration frequency for high quality of measurements
- Non-volatile memory holds up to 1000 data points with date and time in automatic (Data logger) o manual and send to PC or printer
- Stirrer plate included with direct control from pH 80 with independent keys integrated on case or removable with 50cm long cable
- 2 x USB ports to transfer data at PC or Printer and for external keyboard
- Software PC-Link to download data (included)
- Incl. Buffer solutions and accessories
- IP 54 rated housing handles splashes

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Multifunctional bench instrument for measuring pH/ mV, REDOX, Conductivity, TDS and temperature inclusive Polymer pH electrode, conductivity cell 2301T with ATC, magnetical stirrer and accessories.

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- Display pH/mV and Temperature, COND/ TDS/Salinity and Temperature simultaneously
- Selectable Resolution 0,1 / 0,01 pH
- Auto range for COND
- Calibration for up to 3 pH points with autom. recognition of buffer solution (USA / NIST)
- manual calibration for up to 2 pH points defined by user
- Calibration for up to 4 COND points with automatic recognition of buffer solution
- Manual calibration for 1 COND point defined by user
- TDS measuring with selectable factor
- Large 2-color display with indication of Calibration Buffers and measurement's stability
- Stirrer plate included with direct control from PC 80 with independent keys integrated on case or removable with 50cm long cable
- Incl. Buffer solutions and accessories
- IP 54 rated housing handles splashes

Product features

- Display pH/mV and Temperature, COND/ TDS/Salinity and Temperature simultaneously
- Selectable Resolution 0,1 / 0,01 pH
- Auto range for COND
- Calibration for up to 3 pH points with autom. recognition of buffer solution (USA / NIST)
- manual calibration for up to 2 pH points defined by user
- Calibration for up to 4 COND points with automatic recognition of buffer solution
- Manual calibration for 1 COND point defined by user
- TDS measuring with selectable factor
- Large 2-color display with indication of Calibration Buffers and measurement's stability
- Stirrer plate included with direct control from PC 80 with independent keys integrated on case or removable with 50cm long cable
- Incl. Buffer solutions and accessories
- IP 54 rated housing handles splashes

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Multifunctional bench instrument for measuring pH/ mV, REDOX, Conductivity, TDS and temperature inclusive magnetical stirrer and accessories..

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5040-0224

Multifunctional bench instrument for measuring pH/ mV, REDOX, Conductivity, TDS, Resistivity and temperature inclusive STANDARD pH electrode, conductivity cell glass VPT 80/1 with ATC, magnetical stirrer and accessories

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Product features

- Display pH/mV and Temperature, COND/ TDS/Salinity/Resistivity and Temperature simultaneously
- Selectable Resolution 0,1 / 0,01 / 0,001 pH
- Auto range for COND
- Calibration for up to 5 pH points with automatic recognition of buffer solution (USA / NIST / DIN)
- manual calibration for up to 5 pH points defined by user
- Calibration for up to 4 COND points with automatic recognition of buffer solution
- Manual calibration for 1 COND point defined by user
- TDS measuring with selectable factor
- GLP and IQ-OQ manual available
- Large color display with indication of Calibration Buffers and measurement's stability
- Memory with re-call of last calibration date and time and information of electrode's life
- CAL DUE: set of calibration frequency for high quality of measurements
- Non-volatile memory holds up to 1000 data points with date and time in automatic (Data logger) o manual and send to PC or printer
- 2 x USB ports to transfer data at PC or Printer and for external keyboard
- Software PC-Link to download data (included)
- Stirrer plate included with direct control from PC 80 with independent keys integrated on case or removable with 50cm long cable
- Incl. Buffer solutions and accessories
- IP 54 rated housing handles splashes

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5040-0227

Multifunctional bench instrument for measuring pH/ mV, REDOX, Conductivity, TDS, Resistivity and temperature inclusive magnetical stirrer and accessories

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Product features

- Stir speed 0...2500 rpm
- Stir capacity up to 1000ml
- Working temperature +5...+40°C
- Maximum stir bar length up to 30mm
- 220 V Power supply
- Incl. rechargeable lithium battery
- Weight 150 g
- IP 54

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Areas of application

Perfect for mobile purposes in laboratory or in field, specially in terms of the implemented lithium accu

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XS electrodes Standard and Flow are provided with a quality certificate attesting and documents the final inspection carried out on each single electrode. Also on the certificate and label shows the serial number of the electrode for a correct traceability (GLP) of the relations of the calibration of the measurement chain.

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Elektrodenverpackung / Electrode box

interne Seriennummer / Indelible internal serial number

XS Sensorelektroden / Sensor electrodes

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Very famous Swiss producer has developed a complete line laboratory pH, Redox and conductivity electrodes, with glass and plastic body. The new XS sensor are available with S7 head or cable fix with BNC or DIN connector plug. For food application many models are equipped with temperature sensor inside for fast temperature compensation.

Some models are designed to help in daily work in GLP condition with indelible internal serial number and certificate are provided with quality test. Special polymer electrolyte combined with high performance reference system generate a complete range of electrodes for food application in special part of this catalogue.

Blue inner buffer provides clear visual of contact with internal pH glass membrane. All XS sensors line are connectable at all common pHmeters brand in the market. Our customer support are available for every special application request or support for right type electrode indication.

XS sensor for pH is not only pH electrode is perfect combination for pH measurement with XS pH meters and XS buffer solutions.

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
XS Leitfähigkeitsmesszellen / XS Conductivity cells

Modell / model	Bereich COND / Range COND	Temperatur °C / Temperature °C	Anschluss / Connection	Hinweis / Notes	Artikel Nr. / Ord.No.	€
VP 8	10 µS... 10 mS	0°C...+80 °C	Festkabel / BNC / Fix cable / BNC	Zellenkonstante C = 1 / cell constant C = 1	6030-0080 ③	114,00
VPT 80/1	10 µS... 10 mS	0°C...+80 °C	Festkabel / BNC / Fix cable / BNC	Zellenkonstante C = 1 Mit Sensor NTC 30KΩ Stecker RCA Cinch / cell constant C = 1 With sensor NTC 30KΩ plug RCA Cinch	6030-0081 ③	148,00

Material / material	Glas / Glass	
Temperaturfühler / Temperature probe	Nein / Ja in VPT 80/1 / No/Yes in VPT 80/1	
Diaphragma / diaphragm	Platin / Platinum	6030-0080 6030-0081
Länge / length	120 mm	
Durchmesser / diameter	12 mm	

Modell / model	Bereich COND / Range COND	Temperatur °C / Temperature °C	Anschluss / Connection	Hinweis / Notes	Artikel Nr. / Ord.No.	€
VPT 51/01	0,1 µS... 1 mS	0°C...+60 °C	Festkabel / BNC / Fix cable / BNC	Zellenkonstante C = 0,1 mit Sensor NTC 30KΩ Stecker RCA Cinch Reinstwasser / cell constant C = 0,1 With sensor NTC 30KΩ plug RCA Cinch ultrapure water	6030-0082 ③	227,00

Material / material	PVC / Stahl / PVC / Steel	
Temperaturfühler / Temperature probe	JA / Yes	
Diaphragma / diaphragm	Stahl / Steel	6030-0082
Länge / length	120 mm	
Durchmesser / diameter	12 mm	

XS Leitfähigkeitsmesszellen / XS Conductivity cells

Modell / model	Bereich COND / Range COND	Temperatur °C / Temperature °C	Anschluss / Connection	Hinweis / Notes	Artikel Nr. / Ord.No.	€
VPT 80/10	100 µS... 200 mS	0°C...+80 °C	Festkabel / BNC / Fix cable / BNC	Zellenkonstante C = 10 Mit Sensor NTC 30KΩ Stecker RCA Cinch für Meerwasser geeignet/ cell constant C=10 With sensor NTC 30KΩ plug RCA Cinch sea water	6030-0083 ③	221,00

Material / material	Glas / Glass
Temperaturfühler / Temperature probe	JA / Yes
Diaphragma / diaphragm	Platin / Platinum
Länge / length	120 mm
Durchmesser / diameter	12 mm

Modell / model	Bereich COND / Range COND	Temperatur °C / Temperature °C	Anschluss / Connection	Hinweis / Notes	Artikel Nr. / Ord.No.	€
VPT 80/4	10 µS... 200 mS	0°C...+60 °C	Festkabel / multi Stift / Fix cable / multipin	Zellenkonstante C = 0,5 Mit Sensor NTC 30KΩ Stecker Multistift / cell constant C = 0,5 With sensor NTC 30KΩ plug multipin	6030-0084 ③	242,00

Material / material	PVC
Temperaturfühler / Temperature probe	JA / Yes
Diaphragma / diaphragm	4 Graphitplatte / 4 graphite plate
Länge / length	120 mm
Durchmesser / diameter	12 mm

Modell / model	Bereich COND / Range COND	Temperatur °C / Temperature °C	Anschluss / Connection	Hinweis / Notes	Artikel Nr. / Ord.No.	€
2301T	10 µS... 200 mS	0°C...+60 °C	Festkabel / BNC / Fix cable / BNC	Zellenkonstante C = 1 Mit Sensor NTC 30KΩ Stecker RCA Cinch / cell constant C = 1 With sensor NTC 30KΩ plug RCA Cinch	6030-0085 ③	128,00

Material / material	PVC / Stahl / PVC / Steel
Temperaturfühler / Temperature probe	JA / Yes
Diaphragma / diaphragm	Stahl / Steel
Länge / length	120 mm
Durchmesser / diameter	12 mm

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Pufferlösungen in Patent-Flasche / buffer solution N.I.S.T- Certificate in patent bottle

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“The calibration of pH measuring system is one of the most important factors. To ensure a high measuring accuracy one has to follow the correct calibration procedure, use new and accurate buffer solutions and wait long enough for stable indication values of the electrode assembly to be calibrated. No pH measurement can be more accurate than its calibration accuracy” (International of basic and General Terms in Metrology)
If you have a good electrode, but you don't calibrate this with buffer solutions with good metrological traceability, You can't obtain reliable results. The XS buffer solutions come to celebrate the Third Perfect Combination.

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1. General Facts

- 1.1. (Conflicting business conditions, written form). We only conclude to the following terms, which shall apply to our other business transactions with the customer without further notice. Conflicting terms are not accepted even without explicit objection. Agreements with us, delivery times, guaranteed qualities and warranty or obligation costs are only valid if we issue a written confirmation to the customer. The following terms and conditions can also be found in the remarks printed in the appendix of the price list.
- 1.2. (Offers, subject of change, data storage). Our offers are subject to change. Technical improvements to our products can be performed without notice. We may store and process relevant contract data on our processing systems.
- 1.3. (Offsetting, withholding). Except of undisputable counterclaims, offsetting or withholding by the customer is not permitted.
- 1.4. (rush orders/small orders). The order confirmation is the same as the invoice confirmation (attached carbon copy) for delivery within 8 days or order value up to 500 Euro.
- 1.5. (place of execution, jurisdiction, governing law). Place of execution is our plant site in Wertheim Reicholzheim, jurisdiction over registered traders is Wertheim / Mosbach. The German law is applicable to the exclusion of the United Nations Convention on Contracts for the International Sale of Goods (CISG).

2. Risk, shipping costs, number of items, recall, non-acceptance

- 2.1. The risk is transferred onto the customer the moment the products to-be-delivered leave our plant. The customer carries transport and insurance costs until the place of delivery as well as the costs for specialized packaging. A specialized package is every non-standard package exceeding the standard packaging we offer to the customers.
- 2.2. Excess and short delivery of +/- 10% are allowed for agreed quantities. For call orders, the total amount must be collected within 12 months, whereby the rest of the released quantity can be shipped to the customer after the deadline.
- 2.3. If the customer does not accept goods ready-for dispatch or already shipped products in due time, we have the right to place the goods in storage while maintaining our claim for performance or to sell them according to 5.5.

3. Delivery, delay

- 3.1. Delivery time is only an approximate estimation, counting from the moment the customer receives our order confirmation, whereby delivery time starts earliest after clarifying technical questions and collecting installments as well related documents from the customer and ends with the shipment. Exceeding the deadline up to two weeks does not have legal consequences.
- 3.2. Acts of God and situations beyond our responsibility like strikes, lockouts, breakdowns, raw material and equipment shortages and delayed delivery by suppliers or customer-requested additional or modified services, prolong the delivery time accordingly and release us from our obligation to deliver if delivery is impossible.
- 3.3. At fixed time arrangements we can only default through a reminder of the customer. In order to make claims resulting from the delay, the customer needs to set a reasonable grace period after occurrence of default. We are only liable for delay damages caused by us or our directors, employees and auxiliary persons intentionally or through gross negligence. Liability is limited to the foreseeable damages at contract conclusion.

4. Payment terms, price changes, return compensation

- 4.1. Prices are quoted ex works. Invoices are payable within one month after sending net or in 10 days with a 2% discount. We only take drafts and checks at the expense of the customer on account of performance. If there is doubt as to the creditworthiness of the customer, we can demand payment in cash, cash in advance or security deposits.
- 4.2. In accordance with § 315 BGB, we may demand a surcharge in proportion to our increase in costs (including tax increases) in reasonable discretion between contract conclusion and delivery. The customer may cancel in case of price increases of more than 15%. Call deliveries are priced according to our daily price. Changes of agreed dimensions, numbers, etc. can result in recalculation.
- 4.3. A compensation of 15% of the invoiced value is payable if it is agreed that goods free-of-defects are returned.

5. Retention of title, assignment in advance

- 5.1. The goods to-be-delivered remain our property until their full payment without restrictions. The retention of title also applies to the redemption of still outstanding claims with the customer if he is a registered trader. The customer may combine the goods to-be-delivered with other items before payment only if they are not burdened with rights of third parties.
- 5.2. Modifying the goods to-be-delivered is free of charge for us. If we lose our ownership of the delivered goods with installation, we will become co-owners of the newly-built products in proportion to the value of the goods to-be-delivered. The customer stores conditional goods for us at no charge.
- 5.3. The customer may sell our conditional goods (paragraph 5.1 and 5.2) in the ordinary course of business only if his claims from the sale have not previously been assigned, seized or charged otherwise or can be set-off with counter-claims of his

customer and if no late payment to us has occurred. In case of conflicting assignments in advance in accordance with our scope of supply, he assigns all claims to us arising from the sale to his customers or from advance payment of such claims to factoring-banks in the amount of our claims (see paragraph 5.1) as a security deposit. In factoring, the customer may sell conditional goods only if the factor knows about the assignment to us and pays the assigned partial value (paid to the exclusion of a collection authority of the customer) directly to us.

5.4. The customer must keep the incoming payments for the assignment separate for us and use them to settle our claims. We float claims of more than 20% of excess collateral on request.

5.5. In case of late payment, including payment of goods from previous deliveries, we can demand that the customer returns all remaining conditional goods and removes our possibly installed products, also we can forbid the resale of these products; furthermore we can request the disclosure of the assigned claims according to paragraph 5.3 and forbid its collection. Redemption or recovery of conditional goods takes place with respect to the estimate value and appropriate credit note.

6. Warranty, compensation for damages, spare parts

- 6.1. Only the warranted characteristics and other commitments we expressed towards the customer in written form are binding. Information in brochures and instruction manuals or references to industrial norms do not justify warranted characteristics or acceptance of specific obligations to meet claims. If the customer needs the goods for special purposes, he must check its suitability for the intended usage - also in terms of product safety - in advance, especially if the products meet all relevant technical or regulatory requirements. Without prior testing, claims for compensation due to damages resulting from unsuitability are invalid. In case of material or design requirements of the customer we are not liable for the suitability or eligibility of the desired materials or designs and thus have no particular duty to review. For non-calibratable products, variations within the double-error limits are still in accordance with the contract.
- 6.2. The customer loses his warranty and compensation claims resulting from patent defects or absence of guaranteed properties if he does not check the delivered goods for product safety immediately after receipt, at the latest before processing, consumption, use, installation or resale and hands in his complaint within 10 days in writing. In addition such claims (even with hidden defects) become time-barred within six months after delivery.
- 6.3. In case of legitimate complaints and after the customer has set a reasonable period, we are according to paragraph 6.2. only required to repair, to swap or to additionally supply the delivered goods or removable product parts free of charge according to our own judgment if the products are defective due to factors before risk transfer such as faulty production, low-quality material or poor workmanship. Only in case of unfounded refusal, failure or inability to ensure the above measures can the customer demand a change or a reduction and in the absence of guaranteed properties compensation for damages. We are only liable for consequential damages if the customer explicitly states the possible dangers of these damages at contract conclusion and we accept a particular cost obligation on this risk.
- 6.4. Other claims of the customer such as a positive breach, tort (including product liability) or other legal grounds (such as advice, instruction manuals, maintenance, fault at contract conclusion) count against us only if deliberate intent or gross negligence can be found in our behaviour and that of our directors or auxiliary persons. The claims are restricted to the extent being stated by the customer and being predictable by us at contract conclusion.
- 6.5. Warranty or indemnity claims are excluded if they are attributable to improper handling, maintenance, operation or treatment by the customer or a third party or due to normal wear and tear (especially with wearing parts) or shipping damage.
- 6.6. Warranty and claims for replacement items and other remedial actions are also regulated according to these terms and become time-barred with the end of the warranty period for the original object.
- 6.7. If we have an obligation to maintain spare parts, this is limited for a period of five years after delivery.

7. Industrial property rights, tools, secrecy

- 7.1. For forms, patterns, designs, technical documents, quotations or offers provided by us, we reserve the property and all intellectual property rights and copyrights. The customer may use it only in the agreed manner. Without our written agreement, the customer is not allowed to produce the goods referred to in the contract himself or with the help of others.
- 7.2. If we deliver products according to the buyer's drawings, models and patterns, the customer is responsible that no industrial property rights and other rights of third parties are violated through their production and delivery; furthermore the customer must compensate us for all damages resulting from such abuses.
- 7.3. Forms, tools or other devices produced or provided by us shall remain our property even if the customer has taken over their costs partially or completely.
- 7.4. The customer must keep all non-public knowledge from the business relationship with us secret to third parties.

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We are standing for competence and know-how when it comes to measuring temperature, humidity, air flow and pressure value. Our measuring instruments are manufactured on modern production lines and are carefully assembled. In our well-equipped laboratory, the instruments are matched to comparative references that are traceable back to national standards defined by the German PTB or similar European authorities. Our measuring probes are of high quality, subject to constant testing and accurate controls. This is a guarantee for highly accurate and reliable products. Helping customers with professional advice is part of our company identity. Please feel free to ask for more information!

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DOSTMANN electronic GmbH